3/97-A-Sf-


Fig. 1

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Feldemissions-Mikroskop Field Emission Microscope

1 Safety notes

- Please follow the instructions given in Section 4 "Operation" (pages 4 to 6) to the letter so that faulty operation, which can cause damage to the instrument, can be prevented.
- Use only the devices recommended for this experiment set-up.
This applies in particular for the high-voltage supply and the high-voltage proof heating-voltage supply (see section 4.1).
- The field emission microscope is an interference radiator within the meaning of German x-ray regulations.
Operating the apparatus with the specified high-voltage power supply units 10 kV (521 70 and 255 37*) ensures that the dose rate of the emitted x-ray radiation remains within the permissible maximum value. This applies for a field emission microscope with undamaged tungsten tip, which produces an image of the type shown in Fig. 4/5.
A field emission microscope with a defective tungsten tip (no uniform image as shown in Fig. 4/5, but sharply defined luminance phenomena on the screen) must not be operated with anode voltages greater than 5 kV. Under unfavorable conditions, excessive dose rates may result due to the x-ray radiation generated in the field emission microscope.
- Always clean the cathode during experiment intervals by heating it in the prescribed manner for a period of 60 s.
Severe contamination of the cathode caused by absorbed residual gases will shorten the life of the tube.
- **Barium heating:**
7.5 A to 8 A (until a dull red glow appears)
Place the ammeter on an insulated surface, e.g. glass plate.
Intensive barium vaporization makes the tube unusable.
- **Cathode heating:**
1.6 A to 1.8 A; maximum 1.9 A
Use a class 1.5 (or better) meter; measuring range ≤ 3 A
For maximum cathode life, keep the heating current as low as possible.
Currents over 1.9 A can vaporize the tungsten. This destroys the cathode immediately.
- Avoid sudden applications of high voltage. Make sure that the sockets make good contact. Protect the fluorescent screen from direct sunlight.
- Please bear in mind that the tube is evacuated; mechanical damage to the tube could be hazardous due to flying glass splinters.

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The field emission electron microscope renders the surface of a tungsten monocrystal visible at a scale of magnification of approx. 500,000, by forming an image of the point-to-point differences of the electron emissivity of the surface. With this very high magnification and a resolving power of 2 nm to 3 nm, it is possible to observe individual barium atoms deposited on the tungsten surface, and their thermal agitation.

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Fig. 2

Schematic diagram of the field emission microscope:

- ① Cathode (tungsten tip)
- ② Anode with barium supply in ring-shaped trough
- ③ Heating for barium supply
- ④ fluorescent screen

Fig. 2 shows a field emission electron microscope in schematic representation. A lightbulb-like glass vessel contains at the center of the internally mounted fluorescent screen ④ the cathode ① shaped as a fine point. In Fig. 2 this tip is shown enlarged. It forms one end of a tungsten wire bent into a loop, which is connected to the lamp socket (E 27). The anode ② with the ring-shaped trough for the barium (Ba) is situated laterally to the tip. The anode and barium heating are connected to a lamp socket (E 14). The glass bulb is evacuated to produce an extremely high vacuum, and fused.

The tip of the tungsten cathode has a diameter of less than one wavelength of visible light. Such fine tips are made by etching processes, whereby their initially irregular shape is rounded off by baking them in an ultra-high vacuum. After this treatment, the monocrystal tip no longer appears as a sharp-edged polyhedron, but as a rounded-off form, in agreement with theory. The flattened regions here represent the so-called equilibrium shape faces: (011), (001), (112). Due to the flatness and the relatively high electron work function these regions appear as dark patches in the microscope image of the clean cathode tip.

Technical data

Magnification:	approx. 500,000x
Resolving power:	2 nm to 3 nm
Radius of curvature of the tungsten tip:	0.1 μm to 0.2 μm
Dia. of bulb:	10 cm
Pressure:	approx. 10^{-7} hPa ($\cong 10^{-7}$ mbar)
Anode voltage U_A :	typically 4 kV to 8 kV
Anode current I_A :	approx. 10 μA
Current I_W for heating the tungsten tip:	1.6 A to 1.8 A; max. 1.9 A
Current I_B for heating the barium:	7.5 A to 8 A max.

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3 How the microscope works

The following two points constitute the essential difference between a field emission microscope and an electron microscope working with electrostatic or electromagnetic lenses. The electron source in the field emission microscope is not an incandescent wire, but an extraordinarily fine rounded tungsten tip.

- What is imaged is not an object brought into the path of the beam, but the electrons themselves emitted from the tip surface, by virtue of the field about this tip.
- Between the tungsten tip and the concentrically arranged glass section of the tube, a field with spherical symmetry is established, whose lines of force proceed radially outwards. Near the cathode surface the intensity of the field attains very high values.

Thus, in the case of a cathode radius of $R = 10^{-6}$ m and a voltage $U = 10^4$ V, the field strength would amount to 10^{10} V/m. However, such field strengths cannot be maintained under normal conditions. In air at atmospheric pressure, the breakdown field strength of air, which is approx. $2 \cdot 10^6$ V/m to $3 \cdot 10^6$ V/m, sets a limit to the strength of the field at the surface of the sphere. Even at reduced pressure this strength cannot be exceeded as long as a spontaneous discharge is possible with the residual gases present. Only in a tube evacuated to 10^{-7} bar by modern means, in which spontaneous discharges no longer occur, can the electric field strength around a fine tip with hemispherical form be increased to high values. If the intensity at the negatively poled tip becomes of the order of some 10^9 V/m, conduction electrons can escape from the metal into the vacuum. This emergence of electrons from the metal surface is a physical process that can only be explained by the concepts of wave mechanics (tunnel effect).

Such electrons issuing from the fine cathode tip by field emission move toward the lines of force of the electric field. In the proximity of the tip, these proceed outwards approximately radially from the hemispherical cap. In the extremely high vacuum, the electrons travel to the fluorescent screen in a radial direction practically without undergoing any collisions on their way.

Two electrons leaving the hemispherical tip of the cathode at a linear distance Θ will impinge on the fluorescent screen at a linear distance Δ . The magnification ratio $G = \Delta/\Theta$ of the field emission microscope is simply given by the quotient of the radii of the glass bulb and tip:

$$G = \frac{\Delta}{\Theta} = \frac{R}{r}$$

If the radius of the hemisphere r of the tip is 0.1 to 0.2 μm , which is possible in practical terms, and the radius of the spherical part of the bulb is $R = 5$ cm, then $G = 5 \cdot 10^5$ to $2.5 \cdot 10^5$. The largest magnification of which the field emission microscope is capable is approx. 1,000,000x. This substantially exceeds the magnification possible with a conventional electron microscope.

Under these conditions, distances on the fluorescent screen of the field emission microscope of about 0.5 mm correspond to separate points on the hemispherical tip that have a distance of only $1 \cdot 10^{-9}$ m to $2 \cdot 10^{-9}$ m. But such distances are already in the order of magnitude of atomic dimensions. Thus, e.g., the spherical diameter of a barium atom is approx. $0.4 \cdot 10^{-9}$ m. Therefore, the field emission microscope can reproduce such atomic dimensions as images on the fluorescent screen.

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4 Operation

Important: Please observe the safety notes provided on page 1!

4.1 Accessories required; circuit (see Fig. 3)

High voltage:	
1 High-voltage power supply unit, 10 kV	521 70
or	
1 High-voltage power supply unit, 10 kV	522 37*
Heating voltages:	
1 high-voltage proof, variable extra-low voltage source for barium heating, loading capacity up to 10 A;	
1 U-core with yoke	562 11
1 Clamping device	562 12
1 Extra-low voltage coil	562 18
at 230 V mains voltage:	
1 Coil with 1000 turns	562 15
1 Rheostat, 1000 Ω	537 22
for 115 V mains voltage:	
1 Coil with 500 turns	562 14
1 Rheostat, 320 Ω	537 23
1 Variable extra low-voltage source, loading capacity up to 2 A DC, for cathode heating, e.g. Variable extra-low voltage transformer S	
	521 35
2 Ammeters, measuring range 3 A DC and 10 A DC, class 1.5, e.g. Demo-multimeter	
	531 911
1 Glass plate (insulating surface for ammeter of I _{BA}) e.g.	
	664 196
1 Safety connecting box	
	502 06
1 E 27 socket on rod	
	451 18
1 E 14 socket	
	505 68
1 Stand base, length of sides 20 cm	
	300 02

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4.2 High voltage supply

It is essential that the voltage be applied to the field emission microscope only **after** baking of the tungsten tip is completed and the voltage is **slowly** increased **from zero** until a sufficiently bright image is achieved in the dark room.

Each time before baking out the cathode, set the high voltage back to zero and wait until the output voltage drops to zero (max. 30 seconds for 10 kV power supply 521 70). Ground the

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output if necessary.

4.3 Heating voltage supply

The high-voltage supply for barium heating must have a sufficiently high insulation resistance. The Leybold demonstration transformer is very suitable whereas other transformers generally are not. Ammeters to observe the barium heating current I_{BA} must be arranged on a well-insulated surface, e.g. on a glass plate (664 196).

4.3.1 Baking out the cathode

The cathode must be baked out **when no high voltage is applied to the tube ($U_A = 0$)**. The purpose of baking out is to clean the tungsten tip of all foreign atoms absorbed.

The cathode must be baked out:

- before experimenting with the tube for the first time on any particular day or
- after the end of a demonstration with the barium-coated cathode,
- and whenever this appears to be necessary.

The bake-out operation requires a current of 1.6 to 1.8 A which **must be increased continuously from zero up to the maximum value** and should be maintained at full strength for approx. 1 min.

As the very sensitive tip is easily destroyed, it must be baked out with the utmost care with the following points taken into particular account:

Never apply high voltage and full cathode heating current simultaneously. Eliminate residual charge of the source of high voltage prior to each bake-out.

Never apply the heating current in full strength, but increase gradually under constant observation of the ammeter for I_W ; never exceed the maximum value of 1.9 A.

In the case of abnormal luminous phenomena, switch off high voltage at once and bake out cathode.

If, in the course of operation or when starting up, the luminous area on the screen contracts into a small, bright spot, then the cathode tip has not been baked out long enough or intensely enough prior to the experiment. That is why the cathode must be baked out in the above manner after switching off high voltage. Should this disturbance recur persistently, the cathode heating current may be applied for up to 5 minutes.

Electrical conduction between the anode terminal and the fluorescent screen (acting as the anode) is achieved by vacuum deposition of the barium film within the glass bulb, and via the secondary electrons. When the apparatus has not been put into operation for a longer period of time, it may first be necessary after bake-out of the tungsten tip, to vaporize some barium in order to obtain an image on the screen.

4.3.2 Vacuum deposition of barium

The heating current I_{Ba} is set to approx. 7.5 to 8.0 A starting from 0 (adjust, if necessary). The heating time is 30 s to 60 s.

During the heating procedure observe the ring-shaped trough with barium and the image on the screen. After about 30 s the dark-red glow of the ring becomes visible in the section at the side of the glass bulb, and shortly afterwards a slight movement within the image on the green fluorescent screen. After an additional few seconds many bright points light up, clustered around the large dark points of the lattice structure. The barium heating should now be switched off immediately.


The barium reserve lasts for more than 50 vaporization cycles. After barium has been deposited on the tip, the cathode can be slightly heated while high voltage is on, to observe the thermal agitation of the barium atoms on it. However, the heating voltage should only be increased sufficiently to set the barium atoms in vigorous flickering, the **cathode must never show more than a dull red glow.**

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5 Experiments

5.1 Image of the tungsten tip

The fine tungsten tip of only 0.1 to 0.2 μm radius consists of a monocrystal. This metal has a cubic body-centered lattice, the unit cell of which is represented in Fig. 6. The edge of the unit cell is 0.316 nm long. On account of the manner in which the crystal was prepared, the alignment of the crystallographic axes is such that one diagonal of the unit cell lies in the longitudinal direction of the tip and intersects the zenith of the hemisphere. This axis (110) is indicated by a broken line in Fig. 6. The plane perpendicular to it is labelled (110) in accordance with the standard convention for crystallographic indices. As shown in Fig. 4 and 5 only a few electrons emerge from this plane as revealed by the dark patch at the center of the image. In this direction the electron work function is slightly greater. The other, dark patches in the electron image which are arranged symmetrically around the first can also be clearly related to crystal axes and crystal planes. This is represented in a parallel projection in Fig. 7.

The image of the tungsten tip in the field emission microscope opens up an interesting chapter from the crystallography of simple lattices. As follows from the difference between Figs. 4 and 5, even the slightest contamination of the surface causes a noticeable change in the emission conditions of the cathode. These changes can become even more striking when the cathode is heated in the electric field. Then further bright and dark patches appear on the fluorescent screen. However, do not let these interesting phenomena tempt you into heating up the tip to visible glow while high voltage is on, as this can destroy it.

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5.2 Imaging individual atoms

It is particularly important that atoms and molecules of other substances can be deposited on the tungsten surface, and observed there. The residual gas atoms present in the bulb are not suitable for these purposes; the particles must be as large as possible. Barium atoms are particularly suitable. When the high voltage is connected to heat up the barium in the tube to vaporization, some of these atoms will precipitate on the fine, point-shaped cathode; small bright points will then flash up on the familiar image of the tungsten cathode. They settle like snowflakes on the tip in growing quantities. During this experiment, not too much barium should be vaporized; the heating should soon be switched off again. That fact that these luminous points are images of individual barium atoms was proven by Professor Dr. E.W. Müller in numerous experiments and using various methods.

The barium atoms that have been vapor-deposited on the tungsten tip can be seen particularly clearly on the parts that have a darker appearance due to the alignment of the tungsten lattice. Thus they appear, for example, on the central dark patch resulting from the [110]-direction. On the outside, near the boundary of the field of vision, the directions [011], [101], [101] and [011] appear very clearly, their edges being fringed with brightly shining images of barium atoms. But in the bright fields between, the barium atoms can also be clearly recognized.


Fig. 6 Elementarzelle des Wolframgitters
Unit cell of tungsten lattice

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The crystallographic directions of the tungsten lattice on the hemispherical tip of the tungsten cathode of the field emission microscope, seen in parallel projection

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If, now, the current I_W for heating the tungsten cathode is increased very slowly, starting from a value of 0 V, then the motion of the barium atoms on the various crystal faces of the tungsten tip can be observed. The barium atoms at the outer parts of the image on the screen are the first to come into motion, then gradually those in the central regions soon begin altering their positions. However, the lateral motion itself cannot be detected due to the high velocity. All that is seen is the flash as an atom settles on the new site.

If the temperature is increased a bit more, the motion of the barium atoms on the planes characterized by the dark faces becomes so vigorous that it soon becomes impossible to identify individual barium atoms. At this stage, the images of barium atoms that are visible on the bright zones will begin to move.

Due to the vigorous thermal agitation individual images cannot be isolated. A restless, swarming glow is seen. On further increase of the cathode temperature the barium atoms disappear more and more from the bright zones and contract to characteristic spots between the dark patches on the image. At these locations barium atoms will then be seen to swarm in great agitation. A further increase of temperature is not recommended, as, during this last observation, the cathode already exhibits a dull red glow.

The quantity of barium deposited on the cathode can be slightly increased in a further experiment. Here, the high voltage should be switched off during vaporization. The object of this measure is to prevent the electric field from attracting those barium atoms which travel past the tip to the opposite side of the tip.

On switching the high voltage on again a short time after depositing the barium atoms, it will be clearly seen that the barium atoms have only settled on the side facing the barium. Here, several layers of atoms may be deposited one on top of the other. The image on the fluorescent screen will then be characterized by a bright luminous edge. This is due to the monoatomic boundaries of the deposited barium patch.

This larger quantity of barium on the tungsten cathode can be distributed all over the hemispherical tip by means of heating. Again, start with a very low cathode heating voltage and observe the motion of the barium atoms across the uncoated part of the hemispherical tip. During this process, the dark parts of the image of the tungsten cathode generally remain completely dark; this is due to the fact that these sites are not preferred by the barium atoms. On the contrary, the barium tends to occupy the bright parts of the image like a fluctuating granulation, as in the previously described experiment.