SPEKTRAL-LAMPEN

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SPECTRAL LAMPS

For AC operation*) with ballast only!

Operation of the lamp requires a suitable and correctly adjusted ballast. (Lamp Hg/5 is already provided with a built-in resistor). Mainly inductive chokes adjusted to the lamp current are to be used (for DC Sodium Spectral Lamp an ohmic resistor only). To operate different spectral lamps with the same ballast, it is recommended to use a choke provided with an adjustable air gap or with different taps to vary operating current from 1...1,5 amps. The admissible variation of the operating current from the nominal values (see table) is $\pm 5\%$. In case a resistor is used instead of a choke, it must be suitable to absorb at rated current (see table) the difference between supply voltage and operating voltage.

The full radiation output is reached after a **warm-up time** of some minutes (with exception of Neon and Helium Lamps). Restarting of Mercury Spectral Lamps Hg, Hg/1, Hg/2 and Hg/3 demands a cooling time of a few minutes.

Operating position: vertical, base down. (Na/1 also horizontal, Hg/4 and Hg/5 any).

Base: P 28 (medium prefocus). Exceptions: TI/1 B22d (D.C. Med. Bay.), Hg/4 S9/10 and Hg/5 E14 (German Cand.).

Dimensions: overall length (incl. base) max. 155 mm, ($6\frac{3}{32}$ inches); average diameter: 30 mm ($1\frac{3}{16}$ inches); distance light centre to base flange 75 mm ($2\frac{61}{64}$ inches). Exceptions: TI/1 81 mm ($3\frac{3}{16}$ inches) overall length, 20 mm ($1\frac{13}{16}$ inches) diameter, 46 mm ($1\frac{13}{16}$ inches) light centre to top of base pin. Hg/4 52 mm ($2\frac{1}{16}$ inches) overall length, 8 mm ($\frac{5}{16}$ inches) diameter, light centre in lamp centre. Hg/5 115 mm ($4\frac{9}{16}$ inches) overall length, 14 mm ($\frac{9}{16}$ inches) diameter, 74 mm ($2\frac{15}{16}$ inches) light centre to bottom contact of base.

Operating characteristics:

Rated Voltage: 220 V. Minimum supply voltage: 200 V.

Make sure of correct polarity when using DC sodium spectral lamps, anode to be connected to edge of base.

*) with exception of DC Sodium Spectral Lamp.

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