Lenardröhre

Lenard's Tube

Tube de Lenard

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1. Description

This tube is used for observing cathode rays outside the discharge tube. In order to do this, the electrons are allowed to escape through a thin metal foil.

The tube should be evacuated on a gas ballast pump D 2 (110 12 S or 110 13 S). Voltage is obtained from the induction coil (521 11).

The Lenard's tube consists of an elongated glass cylinder with a female ground joint NS 19 L to fit on the pump. It is coated on the outside with black lacquer. At one end is the connection for the cathode, whereas that for the anode is located above the joint. A perforated metal cap is cemented to the front surface opposite to the cathode.

A small fluorescent screen on a shaft and several aluminum foils are supplied with the apparatus.

2. Attachment of the aluminum foil

Before starting the experiment, an aluminum foil must be adhered vacuumtight to the perforated front surface of the tube.

It is commendable to make a paste by dissolving shellac in a little alcohol, and to apply it thinly to the edge of the tube's front surface leaving the

Notes

1. The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
2. The data DC... refer to the experiment descriptions published as a collection of leaflets in the "LEYBOLD PHYSIKALISCHE HANDBLÄTTER" (PHYSICS LEAFLETS).
3. The data and illustrations are not binding in every detail for the design of the apparatus. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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perforations free. Subsequently press on the aluminum foil tightly along the edge making sure that there are no wrinkles in it. Allow the cement to dry slightly and then heat it carefully over a hot-plate or the like until the foil is firmly fixed, without folds, and the alcohol vaporizes. Too much heat should not be applied otherwise the cement between the tube and the metal cylinder will be affected.

As the metal cap is heated rather strongly by the electrons, the drying process can also be carried out on the pump. As described below, properly connect the induction coil, allow the gas discharge to burn for about 15 minutes with the pump running and gas ballast valve open, and tightly press the foil again.

The method of attaching the aluminum foil with viscous Ramsay grease (177 32) does not render a completely reliable sealing, and therefore cannot be recommended.

3. Operation

Operating voltage : approx. 70 kV
 Amperage : max. 1 mA
 Working pressure : 10^{-2} Torr (mm Hg)

The tube should be fitted on the well-greased vacuum connection of the gas ballast pump D 2 by its female ground joint. The cathode and the anode should be connected to the plate and the tip respectively of the correctly poled induction coil (see Directions for Use No. 521 11).

The required vacuum is obtained when the spark flashes over with a distance of 25 to 30 mm between the tip and the plate. The electrons can be demonstrated with the small fluorescent screen up to a distance of 2 cm. The screen must be held with its fluorescent colour facing the aluminum foil.

4. Precautionary measures

Apart from electron rays, X-rays also appear in the Lenard's tube. The dose rate amounts to 1 mr/min. at a distance of 1 m. It is therefore recommended to demonstrate these experiments only for a few minutes at a time, and to keep the pupils at least 2 m away from the tube.

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